

5X Sample Buffer

Cat. No.: MB01015

Description

5X Sample Buffer is used as a tracking dye for SDS PAGE Gels loading. This buffer contains SDS and is suitable for denaturing gel electrophoresis. Dilute 1:4 with sample before loading.

Note: 5X Sample Buffer should be brought to room temperature prior to use. Mix well before loading gel.

General Specifications

Buffer Type:	5×Sample Buffer
Form:	Buffer
Product Size:	5 mL
Gel Compatibility:	ExpressPlus™ PAGE Gels and SurePAGE™ Gels
Concentration:	125mM Tris, 10%SDS, 0.25%BPB, 10% 2-Mercaptoethanol, 50% Glycerol
Shipping Condition:	Room Temperature
Regulatory Statement:	Not for use in diagnostic procedures

Contents and Storage

Contents: 5X Sample Buffer

Storage: -20°C for 24 months

(주) 바이오디

(우.14322) 경기도 광명시 하안로 60 광명SK테크노파크 A동 1102호
T_ 02.6264.3399 F_ 02.6264.3400 E_ order@bio-d.co.kr

www.bio-d.co.kr