

Transfer Buffer Powder

Description

Transfer Buffer Powder is used for ExpressPlus™ PAGE Gels and SurePAGE™ Gels transfer. Using proprietary techniques, Transfer Buffer Powder is made to have long shelf life, high resolution, smaller volume.

Transfer Buffer Powder is optimized for the transfer of proteins from ExpressPlus™ PAGE Gels and SurePAGE™ Gels to membranes for western blotting. When combined with 10% methanol or ethanol, Transfer Buffer can be used for wet transfer.

General Specifications

Buffer Type:	Transfer Buffer
Form	Powder
Formulation:	25 mM Tris Base, 25mM Bicine
Concentrated	1X
Product Size:	1Box (10/PK)
Gel Compatibility:	ExpressPlus™ PAGE Gels and SurePAGE™ Gels
Shipping Condition:	Room Temperature
Regulatory Statement:	For Research Use Only. Not for use in diagnostic procedures.

Instructions: Reconstitute with 900 ml H₂O, then add 100 ml of methanol or ethanol per bag.

Contents and Storage

Contents: Transfer Buffer Powder

Storage: Room Temperature for 12 months

(주) 바이오디

(우.14322) 경기도 광명시 하안로 60 광명SK테크노파크 A동 1102호
T_ 02.6264.3399 F_ 02.6264.3400 E_ order@bio-d.co.kr

 검색창 “바이오디” 검색
www.bio-d.co.kr

 **bio** *- biosis Design*