

MES SDS Running Buffer Powder

Cat No: M00677

Version 07272016

Product Description

MES SDS Running Buffer is formulated for running ExpressPlus™ and SurePAGE™, Bis-Tris gels. This buffer is ideal for separating small- and medium- sized proteins. For separation of large- to medium- sized proteins, Tris-MOPS SDS buffer (M00138) is recommended.

General Specifications

Buffer Type:	MES SDS Running Buffer
Format:	Powder
Formulation:	50 mM Tris Base, 50 mM MES, 0.1% SDS, 1mM EDTA
Concentration	1X
Product Size:	1 Box (5PK/Box)
Gel Compatibility:	ExpressPlus™ and SurePAGE™, Bis-Tris gels
Shipping Condition:	Room Temperature
Storage:	Room Temperature for 12 months
Regulatory Statement:	For Research Use Only. Not for use in diagnostic procedures.

Instructions

Reconstitute with 1000 ml deionized water to make 1X running buffer per pack of powder.

Do not use acid or base to adjust the pH.

(주) 바이오디

(우.14322) 경기도 광명시 하안로 60 광명SK테크노파크 A동 1102호
T_ 02.6264.3399 F_ 02.6264.3400 E_ order@bio-d.co.kr



검색창 “바이오디” 검색

www.bio-d.co.kr