

Tris-MOPS-SDS Running Buffer Powder

Cat No: M00138

Version 08182016

I	Product Description	1
II	Storage	1

I Product Description

Tris-MOPS-SDS Running Buffer Powder is used for ExpressPlus™ PAGE Gels running. Using proprietary techniques, the powder is made to have long shelf life, smaller volume, high resolution and fast electrophoresis.

Tris-MOPS-SDS Running Buffer Powder is optimized for the running of proteins from ExpressPlus™ PAGE Gels.

Reconstitute with 1000 ml H₂O to make 1X running buffer per bag of powder. Tris-MOPS-SDS Running Buffer can be used for running.

General Specifications

Buffer Type:	Tris-MOPS-SDS Running Buffer
Form	Powder
Formulation:	50 mM Tris base, 50 mM MOPS, 0.1% SDS and 0.1 mM EDTA.
Concentrated	1X
Product Size:	1Box (5/PK)
Gel Compatibility:	ExpressPlus™ PAGE Gels
Shipping Condition:	Room Temperature
Regulatory Statement:	For Research Use Only. Not for use in diagnostic procedures.

II Storage

Room Temperature for 24 months

(주) 바이오디

(우.14322) 경기도 광명시 하안로 60 광명SK테크노파크 A동 1102호
T_ 02.6264.3399 F_ 02.6264.3400 E_ order@bio-d.co.kr

 검색창 “바이오디” 검색
www.bio-d.co.kr

 **bio** *Design*